



NPTX1 Antibody, Biotin conjugated

Product Code	CSB-PA614901LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15818
Immunogen	Recombinant Human Neuronal pentraxin-1 protein (88-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuronal pentraxin-1 (NP1) (Neuronal pentraxin I) (NP-I), NPTX1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPTX1