



CACNA1C Antibody, HRP conjugated

Product Code	CSB-PA615724LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13936
Immunogen	Recombinant Human Voltage-dependent L-type calcium channel subunit alpha-1C protein (1755-1971AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent L-type calcium channel subunit alpha-1C (Calcium channel, L type, alpha-1 polypeptide, isoform 1, cardiac muscle) (Voltage-gated calcium channel subunit alpha Cav1.2), CACNA1C, CACH2 CACN2 CACNL1A1 CCHL1A1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CACNA1C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.