



RIPK1 Antibody, HRP conjugated

Product Code	CSB-PA618785LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13546
Immunogen	Recombinant Human Receptor-interacting serine/threonine-protein kinase 1 protein (328-458AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-interacting serine/threonine-protein kinase 1 (EC 2.7.11.1) (Cell death protein RIP) (Receptor-interacting protein 1) (RIP-1) (Serine/threonine-protein kinase RIP), RIPK1, RIP RIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RIPK1