



RIPK1 Antibody, FITC conjugated

Product Code	CSB-PA618785LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13546
Immunogen	Recombinant Human Receptor-interacting serine/threonine-protein kinase 1 protein (328-458AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-interacting serine/threonine-protein kinase 1 (EC 2.7.11.1) (Cell death protein RIP) (Receptor-interacting protein 1) (RIP-1) (Serine/threonine-protein kinase RIP), RIPK1, RIP RIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RIPK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.