




GAS6 Antibody

Product Code	CSB-PA619872LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14393
Immunogen	Recombinant Human Growth arrest-specific protein 6 protein (70-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth arrest-specific protein 6 (GAS-6) (AXL receptor tyrosine kinase ligand), GAS6, AXLLG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GAS6

Image



Western Blot

Positive WB detected in: 293 whole cell lysate

All lanes: GAS6 antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 80, 75, 69, 45, 42 kDa

Observed band size: 80 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.