



GAS6 Antibody, Biotin conjugated

Product Code	CSB-PA619872LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14393
Immunogen	Recombinant Human Growth arrest-specific protein 6 protein (70-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth arrest-specific protein 6 (GAS-6) (AXL receptor tyrosine kinase ligand), GAS6, AXLLG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GAS6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.