



KISS1 Antibody, FITC conjugated

Product Code	CSB-PA621883LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15726
Immunogen	Recombinant Human Metastasis-suppressor KiSS-1 protein (20-138AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metastasis-suppressor KiSS-1 (Kisspeptin-1) [Cleaved into: Metastin (Kisspeptin-54); Kisspeptin-14; Kisspeptin-13; Kisspeptin-10], KISS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KISS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.