



IL18R1 Antibody, FITC conjugated

Product Code	CSB-PA623807LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13478
Immunogen	Recombinant Human Interleukin-18 receptor 1 protein (22-128AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-18 receptor 1 (IL-18R-1) (IL-18R1) (CD218 antigen-like family member A) (CDw218a) (IL1 receptor-related protein) (IL-1Rrp) (IL1R-rp) (CD antigen CD218a), IL18R1, IL1RRP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL18R1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.