



KMT5B Antibody, FITC conjugated

Product Code	CSB-PA674569LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4FZB7
Immunogen	Recombinant Human Histone-lysine N-methyltransferase KMT5B protein (100-150AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone-lysine N-methyltransferase KMT5B (EC 2.1.1.43) (Lysine N-methyltransferase 5B) (Lysine-specific methyltransferase 5B) (Suppressor of variegation 4-20 homolog 1) (Su(var)4-20 homolog 1) (Suv4-20h1), KMT5B, SUV420H1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KMT5B