



LRRK2 Antibody, Biotin conjugated

Product Code	CSB-PA722493LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5S007
Immunogen	Recombinant Human Leucine-rich repeat serine/threonine-protein kinase 2 protein (2029-2275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat serine/threonine-protein kinase 2 (EC 2.7.11.1) (Dardarin), LRRK2, PARK8
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LRRK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.