



SND1 Antibody, Biotin conjugated

Product Code	CSB-PA745330LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7KZF4
Immunogen	Recombinant Human Staphylococcal nuclease domain-containing protein 1 protein (560-732AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Staphylococcal nuclease domain-containing protein 1 (100 kDa coactivator) (EBNA2 coactivator p100) (Tudor domain-containing protein 11) (p100 co-activator), SND1, TDRD11
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	SND1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.