




CERS6 Antibody

Product Code	CSB-PA751145LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZMG9
Immunogen	Recombinant Human Ceramide synthase 6 protein (324-392AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ceramide synthase 6 (CerS6) (LAG1 longevity assurance homolog 6), CERS6, LASS6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CERS6

Image



Western Blot

Positive WB detected in: HL60 whole cell lysate


All lanes: CERS6 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45, 46 kDa

Observed band size: 45 kDa



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA751145LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.