



SLC16A12 Antibody, Biotin conjugated

Product Code	CSB-PA757957LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZSM3
Immunogen	Recombinant Human Monocarboxylate transporter 12 protein (431-486AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 12 (MCT 12) (Creatine transporter 2) (CRT2) (Solute carrier family 16 member 12), SLC16A12, MCT12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC16A12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.