



MARK2 Antibody, Biotin conjugated

Product Code	CSB-PA762568LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7KZI7
Immunogen	Recombinant Human Serine/threonine-protein kinase MARK2 protein (373-594AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase MARK2 (EC 2.7.11.1) (EC 2.7.11.26) (ELKL motif kinase 1) (EMK-1) (MAP/microtubule affinity-regulating kinase 2) (PAR1 homolog) (PAR1 homolog b) (Par-1b) (Par1b), MARK2, EMK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MARK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.