




TFAP2B Antibody

Product Code	CSB-PA765731LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92481
Immunogen	Recombinant Human Transcription factor AP-2-beta protein (111-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor AP-2-beta (AP2-beta) (Activating enhancer-binding protein 2-beta), TFAP2B
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TFAP2B

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, MCF-7 whole cell lysate

All lanes: TFAP2B antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 51, 52 kDa

Observed band size: 51, 52 kDa