



DTX1 Antibody, FITC conjugated

Product Code	CSB-PA773055LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86Y01
Immunogen	Recombinant Human E3 ubiquitin-protein ligase DTX1 protein (261-410AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase DTX1 (EC 2.3.2.27) (Protein deltex-1) (Deltex1) (hDTX1) (RING-type E3 ubiquitin transferase DTX1), DTX1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DTX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.