



ALYREF Antibody, HRP conjugated

Product Code	CSB-PA774811LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86V81
Immunogen	Recombinant Human THO complex subunit 4 protein (186-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	THO complex subunit 4, Tho4, Ally of AML-1 and LEF-1, Aly/REF export factor, Transcriptional coactivator Aly/REF, bZIP-enhancing factor BEF, ALYREF, ALY, BEF, THOC4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ALYREF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.