



ALYREF Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA774811LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q86V81 |
| Immunogen | Recombinant Human THO complex subunit 4 protein (186-240AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | THO complex subunit 4, Tho4, Ally of AML-1 and LEF-1, Aly/REF export factor, Transcriptional coactivator Aly/REF, bZIP-enhancing factor BEF, ALYREF, ALY, BEF, THOC4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | ALYREF |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |