



SH3TC2 Antibody, Biotin conjugated

Product Code	CSB-PA819482LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TF17
Immunogen	Recombinant Human SH3 domain and tetratricopeptide repeat-containing protein 2 protein (188-427AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain and tetratricopeptide repeat-containing protein 2, SH3TC2, KIAA1985
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SH3TC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.