



OXER1 Antibody, HRP conjugated

Product Code	CSB-PA837442LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDS5
Immunogen	Recombinant Human Oxoeicosanoid receptor 1 protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxoeicosanoid receptor 1 (5-oxo-ETE G-protein coupled receptor) (G-protein coupled receptor 170) (G-protein coupled receptor R527) (G-protein coupled receptor TG1019), OXER1, GPR170 TG1019
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OXER1