



OXER1 Antibody, HRP conjugated

Product Code	CSB-PA837442LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDS5
Immunogen	Recombinant Human Oxoeicosanoid receptor 1 protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxoeicosanoid receptor 1 (5-oxo-ETE G-protein coupled receptor) (G-protein coupled receptor 170) (G-protein coupled receptor R527) (G-protein coupled receptor TG1019), OXER1, GPR170 TG1019
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OXER1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.