



GHSR Antibody, Biotin conjugated

Product Code	CSB-PA842164LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92847
Immunogen	Recombinant Human Growth hormone secretagogue receptor type 1 protein (1-40AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth hormone secretagogue receptor type 1 (GHS-R) (GH-releasing peptide receptor) (GHRP) (Ghrelin receptor), GHSR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GHSR