



NEDD4L Antibody, Biotin conjugated

Product Code	CSB-PA850411LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PU5
Immunogen	Recombinant Human E3 ubiquitin-protein ligase NEDD4-like protein (289-494AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase NEDD4-like (EC 2.3.2.26) (HECT-type E3 ubiquitin transferase NED4L) (NEDD4.2) (Nedd4-2), NEDD4L, KIAA0439 NEDL3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NEDD4L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.