



ACO2 Antibody, FITC conjugated

Product Code	CSB-PA859943LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99798
Immunogen	Recombinant Human Aconitate hydratase, mitochondrial protein (452-676AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aconitate hydratase, mitochondrial (Aconitase) (EC 4.2.1.3) (Citrate hydro-lyase), ACO2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ACO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.