



DSG1 Antibody, HRP conjugated

Product Code	CSB-PA861938LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02413
Immunogen	Recombinant Human Desmoglein-1 protein (384-539AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Desmoglein-1 (Cadherin family member 4) (Desmosomal glycoprotein 1) (DG1) (DGI) (Pemphigus foliaceus antigen), DSG1, CDHF4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DSG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.