



ARHGEF12 Antibody, Biotin conjugated

Product Code	CSB-PA868370LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZN5
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 12 protein (1302-1501AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho guanine nucleotide exchange factor 12 (Leukemia-associated RhoGEF), ARHGEF12, KIAA0382 LARG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGEF12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.