



TUBGCP4 Antibody, FITC conjugated

Product Code	CSB-PA871402LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGJ1
Immunogen	Recombinant Human Gamma-tubulin complex component 4 protein (203-309AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-tubulin complex component 4 (GCP-4) (hGCP4) (Gamma-ring complex protein 76 kDa) (h76p) (hGrip76), TUBGCP4, 76P GCP4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TUBGCP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.