



GOLPH3 Antibody, HRP conjugated

Product Code	CSB-PA872494LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4A6
Immunogen	Recombinant Human Golgi phosphoprotein 3 protein (16-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Golgi phosphoprotein 3 (Coat protein GPP34) (Mitochondrial DNA absence factor) (MIDAS), GOLPH3, GPP34
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GOLPH3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.