



GOLPH3 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA872494LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H4A6 |
| Immunogen | Recombinant Human Golgi phosphoprotein 3 protein (16-150AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Golgi phosphoprotein 3 (Coat protein GPP34) (Mitochondrial DNA absence factor) (MIDAS), GOLPH3, GPP34 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | GOLPH3 |