



CACYBP Antibody, HRP conjugated

Product Code	CSB-PA872536LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB71
Immunogen	Recombinant Human Calcyclin-binding protein (2-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcyclin-binding protein (CacyBP) (hCacyBP) (S100A6-binding protein) (Siah-interacting protein), CACYBP, S100A6BP SIP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CACYBP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.