



JAM3 Antibody, Biotin conjugated

Product Code	CSB-PA874828LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BX67
Immunogen	Recombinant Human Junctional adhesion molecule C protein (119-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Junctional adhesion molecule C (JAM-C) (JAM-2) (Junctional adhesion molecule 3) (JAM-3), JAM3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	JAM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.