




DBF4 Antibody

Product Code	CSB-PA887010LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBU7
Immunogen	Recombinant Human Protein DBF4 homolog A protein (101-281AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein DBF4 homolog A (Activator of S phase kinase) (Chiffon homolog A) (DBF4-type zinc finger-containing protein 1), DBF4, ASK DBF4A ZDBF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DBF4

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA887010LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.