



DBF4 Antibody, HRP conjugated

Product Code	CSB-PA887010LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBU7
Immunogen	Recombinant Human Protein DBF4 homolog A protein (101-281AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein DBF4 homolog A (Activator of S phase kinase) (Chiffon homolog A) (DBF4-type zinc finger-containing protein 1), DBF4, ASK DBF4A ZDBF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DBF4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.