




BCLAF1 Antibody

Product Code	CSB-PA889149LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYF8
Immunogen	Recombinant Human Bcl-2-associated transcription factor 1 protein (218-407AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bcl-2-associated transcription factor 1 (Btf) (BCLAF1 and THRAP3 family member 1), BCLAF1, BTF KIAA0164
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BCLAF1

Image



Immunofluorescent analysis of HeLa cells using CSB-PA889149LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)