




FOXD3 Antibody

Product Code	CSB-PA890729LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJU5
Immunogen	Recombinant Human Forkhead box protein D3 protein (1-140AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein D3 (HNF3/FH transcription factor genesis), FOXD3, HFH2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FOXD3

Image



Western Blot

Positive WB detected in: Rat kidney tissue


All lanes: FOXD3 antibody at 3.2μg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 48 kDa

Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA890729LA01HU at dilution of 1:100



Immunofluorescence staining of HeLa cells with CSB-PA890729LA01HU at 1:133, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).