



NDRG4 Antibody, FITC conjugated

Product Code	CSB-PA890932LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULP0
Immunogen	Recombinant Human Protein NDRG4 protein (210-331AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NDRG4 (Brain development-related molecule 1) (N-myc downstream-regulated gene 4 protein) (Vascular smooth muscle cell-associated protein 8) (SMAP-8), NDRG4, BDM1 KIAA1180
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NDRG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.