




ING1 Antibody

Product Code	CSB-PA891532LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK53
Immunogen	Recombinant Human Inhibitor of growth protein 1 protein (260-353AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibitor of growth protein 1, ING1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ING1

Image



Western Blot

Positive WB detected in: Mouse heart tissue


All lanes: ING1 antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47, 32, 24, 27, 30 kDa

Observed band size: 47 kDa



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA891532LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.