



FZD10 Antibody, Biotin conjugated

Product Code	CSB-PA891782LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULW2
Immunogen	Recombinant Human Frizzled-10 protein (125-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Frizzled-10 (Fz-10) (hFz10) (FzE7) (CD antigen CD350), FZD10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FZD10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.