



FZD10 Antibody, Biotin conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA891782LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9ULW2 |
| Immunogen | Recombinant Human Frizzled-10 protein (125-225AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Frizzled-10 (Fz-10) (hFz10) (FzE7) (CD antigen CD350), FZD10 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | FZD10 |