



LAMP3 Antibody, Biotin conjugated

Product Code	CSB-PA892177LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQV4
Immunogen	Recombinant Human Lysosome-associated membrane glycoprotein 3 protein (28-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysosome-associated membrane glycoprotein 3 (LAMP-3) (Lysosomal-associated membrane protein 3) (DC-lysosome-associated membrane glycoprotein) (DC LAMP) (Protein TSC403) (CD antigen CD208), LAMP3, DCLAMP TSC403
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	LAMP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.