




TAOK2 Antibody

Product Code	CSB-PA892341LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL54
Immunogen	Recombinant Human Serine/threonine-protein kinase TAO2 protein (719-880AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase TAO2 (EC 2.7.11.1) (Kinase from chicken homolog C) (hKFC-C) (Prostate-derived sterile 20-like kinase 1) (PSK-1) (PSK1) (Prostate-derived STE20-like kinase 1) (Thousand and one amino acid protein kinase 2), TAOK2, KIAA0881 MAP3K17 PSK PSK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TAOK2

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA892341LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.