



MTCH2 Antibody, FITC conjugated

Product Code	CSB-PA897586LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6C9
Immunogen	Recombinant Human Mitochondrial carrier homolog 2 protein (31-131AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial carrier homolog 2 (Met-induced mitochondrial protein), MTCH2, MIMP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MTCH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.