



BDKRB2 Antibody, HRP conjugated

Product Code	CSB-PA002652LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30411
Immunogen	Recombinant Human B2 bradykinin receptor protein (336-391AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	B2 bradykinin receptor (B2R) (BK-2 receptor), BDKRB2, BKR2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	BDKRB2