



CA4 Antibody, Biotin conjugated

Product Code	CSB-PA004372LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22748
Immunogen	Recombinant Human Carbonic anhydrase 4 protein (75-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbonic anhydrase 4 (EC 4.2.1.1) (Carbonate dehydratase IV) (Carbonic anhydrase IV) (CA-IV), CA4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.