



# CD3EAP Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA004932LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O15446
<b>Immunogen</b>	Recombinant Human DNA-directed RNA polymerase I subunit RPA34 protein (289-404AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	DNA-directed RNA polymerase I subunit RPA34 (A34.5) (Antisense to ERCC-1 protein) (ASE-1) (CD3-epsilon-associated protein) (CAST) (CD3E-associated protein) (RNA polymerase I-associated factor PAF49), CD3EAP, ASE1 CAST PAF49
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	CD3EAP
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.