



CD3EAP Antibody, Biotin conjugated

Product Code	CSB-PA004932LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15446
Immunogen	Recombinant Human DNA-directed RNA polymerase I subunit RPA34 protein (289-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed RNA polymerase I subunit RPA34 (A34.5) (Antisense to ERCC-1 protein) (ASE-1) (CD3-epsilon-associated protein) (CAST) (CD3E-associated protein) (RNA polymerase I-associated factor PAF49), CD3EAP, ASE1 CAST PAF49
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CD3EAP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.