



CPT2 Antibody, Biotin conjugated

Product Code	CSB-PA005925LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23786
Immunogen	Recombinant Human Carnitine O-palmitoyltransferase 2, mitochondrial protein (388-551AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carnitine O-palmitoyltransferase 2, mitochondrial (EC 2.3.1.21) (Carnitine palmitoyltransferase II) (CPT II), CPT2, CPT1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CPT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.