



CYP4F2 Antibody, Biotin conjugated

Product Code	CSB-PA006451LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78329
Immunogen	Recombinant Human Phylloquinone omega-hydroxylase CYP4F2 protein (268-373AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phylloquinone omega-hydroxylase CYP4F2 (EC 1.14.13.194) (20-hydroxyeicosatetraenoic acid synthase) (20-HETE synthase) (EC 1.14.13.-) (Arachidonic acid omega-hydroxylase) (CYPIVF2) (Cytochrome P450 4F2) (Cytochrome P450-LTB-omega) (Leukotriene-B(4) 20-monooxygenase 1) (Leukotriene-B(4) omega-hydroxylase 1) (EC 1.14.13.30), CYP4F2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CYP4F2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.