



ELK3 Antibody, Biotin conjugated

Product Code	CSB-PA007606LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41970
Immunogen	Recombinant Human ETS domain-containing protein Elk-3 protein (163-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ETS domain-containing protein Elk-3 (ETS-related protein ERP) (ETS-related protein NET) (Serum response factor accessory protein 2) (SAP-2) (SRF accessory protein 2), ELK3, NET SAP2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ELK3