



EOMES Antibody, Biotin conjugated

Product Code	CSB-PA007701LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95936
Immunogen	Recombinant Human Eomesodermin homolog protein (562-676AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eomesodermin homolog (T-box brain protein 2) (T-brain-2) (TBR-2), EOMES, TBR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EOMES